


Instrument Sterilisation Trays


TRANSPAK STACKPAK
MICROPAK SCOPEPAK
MULTIPAK ENDOPAK



WA Distributor:

MICROPAK®

MICROSURGICAL INSTRUMENT TRAYS



- **Protect**
- **Secure**
- **Flexible**
- **Durable**
- **Convenient**

Instrument protection in a complete range of sizes

MicroPak® instrument trays are designed for use with a wide variety of microsurgical and other delicate instruments. These durable, lightweight plastic trays are engineered for maximised protection and ease of handling, with transparent lids for quick instrument identification. Made of rugged, state-of-the-art plastic polymers, the MicroPak is compatible for use with the majority of sterilisation methods currently available.

Cat#	Description	Dimensions		
		Width	Length	Depth
		mm	mm	mm
0100	Base & Lid	40.6	66.0	12.7
0101	Base, Lid & Mat	40.6	66.0	12.7
0111	Base, Lid & Mat	50.8	152.4	19.1
0112	Base, Lid & Mat	63.5	246.4	25.4
0126	Base, Lid & Silicone Bars (1 diamond knife)	50.8	152.4	19.1
0127	Base, Lid & Silicone Bars (2 diamond knife)	50.8	152.4	19.1
0128	Base, Lid & Silicone Bars (3 diamond knife)	76.2	152.4	19.1
0211	Base, Lid & Silicone Mat	76.2	152.4	19.1
1010	Base, Lid	152.4	254	19.1
1011	Base, Lid & Silicone Mat	152.4	254	19.1
1015	Base, Lid & Mat	139.7	241.3	63.5
1020	Base, Lid	152.4	254	38.1
1021	Base, Lid, Insert Tray & 2 Mats	152.4	254	38.1
1023	Base, Lid, Insert Tray & 1 Insert Mat	152.4	254	38.1
1025	Base, Lid & Mat	152.4	254	38.1
2011	Base, Lid & Mat	254	381	19.1
2020	Base & Lid	254	381	38.1
2021	Base, Lid, Insert Tray & 2 Mats	254	381	38.1
2023	Base, Lid, Insert Tray & 1 Insert Mat	254	381	38.1
2024	Base, Lid, Insert Tray & 2 Mats	254	381	38.1
2025	Base, Lid & Mat	254	381	38.1

WA Distributor:

MEDICAL

■ PHONE: 1800 SONITA(766 482)

■ MOBILE: 0423 178 332

■ FAX: 08 9446 2637

SONITA MEDIC

SURGICAL

■ EMAIL: customers@sonitamedic.com.au


■ WEB: www.sonitamedic.com.au

SUPPLIES

■ OFFICE ADDRESS: Suite 76, Ground Floor 50 St. Georges Terrace, Perth WA 6000

SCOPEPAK®

ENDOSCOPE TRAYS



Protects endoscopes during transport, storage and sterilisation.

The ScopePak® offers engineered protection for the handling and storage of costly endoscopes. Available in a variety of sizes to fit a broad range of endoscopes, each tray features silicone holding bars to secure the telescope in place, and a transparent lid for easy identification. Made of rugged state-of-the-art plastic polymers, the ScopePak is engineered to withstand steam sterilization as well as ethylene oxide or cold sterilising solutions.

- Proven Protection - helps prevent damage to endoscopes, minimising expensive repairs
- Secure - endoscopes are fixed in place with unique silicone holding bars
- Flexible - adjustable silicone holding bars accommodate a variety of endoscope sizes
- Durable - engineered with tough, time-tested, highly resilient plastic polymers to withstand the rigors of repeated use and sterilisation
- Convenient - transparent lids allows easy identification of endoscopes
- Versatile- can be used with steam sterilisation, ethylene oxide and cold sterilising solutions

Cat#	Description	Dimensions		
		Width mm	Length mm	Depth mm
9312	ScopePak®, small w/Scope Holder Sets ¹	76.2	281.94	33.02
9317	ScopePak®, medium w/Scope Holder Sets ²	76.2	408.94	33.02
9321	ScopePak®, large w/Scope Holder Sets ²	76.2	510.54	33.02
9326	ScopePak®, extra large w/Scope Holder Sets ³	101.6	660.4	71.12
9331	ScopePak®, extra long w/Scope Holder Sets ³	119.38	769.62	71.12

Notes: 1. Includes 1 ea. #5350 and 1 ea. #5351 silicone scope holder
 2. Includes 2 ea. #5350 and 1 ea. #5351 silicone scope holder
 3. Includes 4 ea. #5353 and 2 ea. #5352 silicone scope holder
 4. Dimensions are inside dimensions.

WA Distributor:

MEDICAL

■ PHONE: 1800 SONITA(766 482)

■ MOBILE: 0423 178 332

■ FAX: 08 9446 2637

SONITA MEDIC

SURGICAL

■ EMAIL: customers@sonitamedic.com.au

■ WEB: www.sonitamedic.com.au

SUPPLIES

■ OFFICE ADDRESS: Suite 76, Ground Floor 50 St. Georges Terrace, Perth WA 6000

MULTIPAK®

INSTRUMENT ORGANISING SYSTEM



Multiples your potential for organisation, sterilisation and savings.

The MultiPak System comprises:

- 4 sizes of Base Units
- 8 sizes of Insert Trays
- 7 sizes of Metal Instrument Racks
- 9 sizes of Lids

You can 'customise' exactly the tray you need.

Multipak® is your unique instrument processing, handling, transport and sterilisation tray system. You can custom-configure the MultiPak to accommodate almost limitless instrument system applications. And the 'a la Carte' method of ordering lets you buy only the pieces you need. MultiPak™ is the most versatile tray system for organising and protecting delicate surgical instrumentation.

The wide array of sizes ensures the proper tray configuration for even the most complicated sets. Multipak trays can be used alone or layered within the base tray in a variety of space saving combinations, like the versatile Insert Tray lids to fit silicone accessories.


This smart design maximises instrument storage while minimising the overall size of the tray, conserving valuable storage and steriliser space.

Organise your sets and reduce repair costs, all in one complete package. Multipak's versatility allows for instrument set modifications with very little, if any, tray modifications. Engineered for use in the most popular modes of sterilisation in use today (Steam, EtO, Plasma, Cold Soaking), that meets your needs and your budget.

Silicone Pin Mats

The Pin Mats provide a convenient method for the positioning and protection of surgical instruments during sterilisation, storage, or transport. The Soft silicone material withstands repeated autoclaving, and mats are available with or without perforations. Mats may be cut to fit existing trays or baskets.

The Fine-pin Mats have closely spaced flexible pins and are recommended for micro surgical instruments. Coarse-pin Mats have medium spaced semi-rigid pins for less delicate instruments.



WA Distributor:

MEDICAL

■ PHONE: 1800 SONITA(766 482)

■ MOBILE: 0423 178 332

■ FAX: 08 9446 2637

SONITA MEDIC


SURGICAL

■ EMAIL: customers@sonitamedic.com.au

■ WEB: www.sonitamedic.com.au

SUPPLIES

■ OFFICE ADDRESS: Suite 76, Ground Floor 50 St. Georges Terrace, Perth WA 6000




These two brochures illustrate all the components available for setting up your trays to your specification:

- Silicone Accessories Catalogue
- System Configuration & User Manual

Cat#	Description	Dimensions*		
		Width	Length	Depth
Base Units and Lids		mm	mm	mm
5000-01-475P	MultiPak Base Unit, Large	233.68	508	116.84
5000-01-375P	MultiPak Base Unit, Large	233.68	508	91.44
5000-01-275P	MultiPak Base Unit, Large	233.68	508	66.04
5000-50-000P	MultiPak Lid, for 50000 Base, Clear	246.38	520.7	43.18
5000-50-175P	MultiPak Elevated Lid, Clear	246.38	520.7	35.56
5006-01-390P	MultiPak Base Unit, Small	139.7	325.12	97.79
5000-02-000P	MultiPak Lid for 5006 Base	147.32	327.66	43.20
Insert Trays and Lids				
5006-03-150P	MultiPak Insert Tray, Small for 5006 Base	129.54	284.48	25.40
5010-03-180P	MultiPak Full Insert Tray, for 5000 Base	220.98	495.3	41.91
5010-50-000P	MultiPak Lid, for 5010 Insert	231.14	505.46	43.18
5012-03-180P	MultiPak Long Insert Tray, for 5000 Base	101.6	495.3	41.91
5012-50-000P	MultiPak Lid, for 5012 Tray	109.22	505.46	43.18
5022-03-180P	MultiPak Insert Tray, for 5000 Base	220.98	241.3	41.91
5022-50-000P	MultiPak Lid, for 5022 Insert Tray, Clear	228.6	231.14	43.18
5122-03-375P	MultiPak Insert Tray, for 5000 Base, (Deep)	220.98	241.3	83.82
5022-50-175P	MultiPak Elevated Lid for 5022, or 5122 Tray	213.36	231.14	35.56
5023-03-180P	MultiPak Insert Tray, for 5000 Base	220.98	322.58	41.91
5023-50-000P	MultiPak Lid, for 5023 Tray, Clear	228.6	330.2	43.18
5123-03-375P	MultiPak Insert Tray, for 5000 Base, (Deep)	220.98	322.58	83.82
5023-50-175P	MultiPak Elevated Lid for 5023, or 5123 Insert Tray	213.36	330.2	35.56
5023-03-E	MultiPak Insert Tray, for 5000 Base	220.98	325.12	41.91
5023-50-E	MultiPak Lid for 5023-50-E, Insert Tray, Clear	223.52	332.74	43.20

* dimensions are for the inside diameter

Part No.	Description	Size mm
0030.04	Coarse Pin Mat, plain	480 x 250 x 18
0032.04	Coarse Pin Mat, perforated	480 x 250 x 18
0031.04	Fine Pin Mat, plain	480 x 250 x 18
0033.04	Fine Pin Mat, perforated	480 x 250 x 18



WA Distributor:

MEDICAL

■ PHONE: 1800 SONITA(766 482)

■ MOBILE: 0423 178 332

■ FAX: 08 9446 2637

SONITA MEDIC

SURGICAL

■ EMAIL: customers@sonitamedic.com.au

■ WEB: www.sonitamedic.com.au

SUPPLIES

■ OFFICE ADDRESS: Suite 76, Ground Floor 50 St. Georges Terrace, Perth WA 6000

STACKPAK®

MICROSURGICAL INSTRUMENT TRAYS





Use shallow heights

Use tall heights

or Combine the two!

Introducing...
A Breakthrough
Tray System That
Really Stacks Up!

Organise your
delicate instrument
in a neat, simple to
use tray system!
Easy as..1..2..3



StackPak® is a revolutionary new way to handle, transport, organise and sterilise your delicate surgical instrumentation. The unique latch design creates a simple way to combine complex instrument sets into one neat package. It makes organisation a snap and you won't outgrow the system. Compatible for use with the majority of today's complex sterilisation methods.

Cat#	Description	Dimensions**		
		Length	Width	Depth
		mm	mm	mm
StackPak Full	Approx. Size 250 x 500mm			
5710-03-180P	Full length tray	495	221	43
5710-03-375P	Deep full tray	495	221	91
5710-50-000P	Full length lid	495	229	5
StackPak 2/3	Approx. Size 250 x 330mm			
5723-03-180P	2/3 length tray	323	221	43
5723-03-375P	Deep 2/3 length tray	323	221	91
5723-50-000P	2/3 length lid	330	229	5
StackPak 1/2	Approx. Size 250 x 228mm			
5722-03-180P	1/2 length tray	241	221	43
5722-03-375P	Deep 1/2 length tray	241	221	91
5722-50-000P	1/2 length lid	249	229	5

** All measurements are inside dimensions

- New Patented Latch design allows you to stack trays while keeping them secure.
- Choose from two standard tray heights designed to accommodate a wide variety of instrument sets (mix and match heights to optimise tray storage capacity).
- Add additional trays if the set expands. Only ONE lid needed to cover it all.
- Modular System is compatible with all existing Multipak accessories.

WA Distributor:

MEDICAL

■ PHONE: 1800 SONITA(766 482)

■ MOBILE: 0423 178 332

■ FAX: 08 9446 2637

SONITA MEDIC

SURGICAL

■ EMAIL: customers@sonitamedic.com.au

■ WEB: www.sonitamedic.com.au


SUPPLIES

■ OFFICE ADDRESS: Suite 76, Ground Floor 50 St. Georges Terrace, Perth WA 6000

ENDOPAK™

INSTRUMENT TRAYS

THE ENDOPAK™ INSTRUMENT TRAY protects and organises your valuable laparoscopy instruments. Custom configurations can be achieved by using the versatile EndoPak™ silicone bar and pin positioning system. The transparent lid makes instruments and sets easy to identify. The EndoPak™ is engineered and manufactured with durable, highly resilient polymers, and the EndoPak™ is made to withstand repeated exposure to steam, temperatures up to 190°C, ethylene oxide, or cold sterilising solutions.



Protection
Convenience
Flexibility
Durability
Versatility
Cost Savings

- Helps to prevent damage to costly and valuable laparoscopy instruments.
- Protects endoscopy instruments during transport, sterilisation, and storage.
- Transparent lid allows for easy identification of instruments and sets.
- Variety of sizes for differing instrument set needs.
- Custom configurations with silicone bar and pin positioning system.
- Engineered with durable, highly resilient state-of-the-art polymers.
- Withstands repeated exposure to steam, temperatures up to 190 °C, ethylene oxide and cold sterilising solutions.
- Helps to minimise unnecessary and expensive instrument repairs and replacement.

ENDPAK™ INSTRUMENT TRAYS & ACCESSORIES

PRODUCT#	DESCRIPTION	INSIDE DIMENSIONS (MM)		
		LENGTH	WIDTH	DEPTH
9230	Tray with lid, pin set & scope holder set	1457	210	44
9230-01	Tray base only for #9230			
9230-02	Lid only for #9230			
9231	Tray with lid, 2 pin sets & scope holder set (not shown)	457	210	89
9231-01	Tray base only for #9231			
9231-02	Lid only for #9231			
9231-03	Second level tray only for #9231			
9240	Tray with lid, pin set & scope holder set	560	210	44
9240-01	Tray base for #9240			
9240-02	Lid only for #9240			
9241	Tray with 2 levels, lid, 2 pin sets & scope holder set	560	210	89
9241-01	Tray base only for #9241			
9241-02	Lid only for #9241			
9241-03	Second level tray only for #9241			
9230-05	Pin set including 16 pins (order 2 for #9231 and #9241)			
9230-06	1 scope holder set			

WA Distributor:

MEDICAL

■ PHONE: 1800 SONITA(766 482)

■ MOBILE: 0423 178 332

■ FAX: 08 9446 2637

SONITA MEDIC

SURGICAL

■ EMAIL: customers@sonitamedic.com.au

■ WEB: www.sonitamedic.com.au

SUPPLIES

■ OFFICE ADDRESS: Suite 76, Ground Floor 50 St. Georges Terrace, Perth WA 6000

TRANSPAK™

INSTRUMENT SOAKING & TRANSPORT TRAY



Now, there's a safer, more secure way to transport contaminated instruments.

The new TransPak™ is engineered to protect personnel from the hazards of contaminated instruments after surgical use.

This strong, lightweight plastic soaking system is designed for minimised spillage and splashing, helping protect hands and eyes from possible infection during transport to cleaning facilities.

Its removable insert tray allows you to transfer contaminated instruments with no hand contact; if necessary, the TransPak™ itself can be transferred directly to a washer/decontaminator.

Transport Contaminated Instruments with Confidence

Safe

- Minimised spillage and splashing protects personnel.
- Removable insert tray allows transfer of contaminated instruments with no hand contact.
- Container clearly identified in Biohazard Red.

Versatile

- Made of special plastic polymers tough enough to be put directly in a washer/decontaminator.
- Use with cold sterilising solutions.
- Engineered to withstand the rigours of repeated use.

Convenient

- Removable insert tray designed to accept MultiPak™ and most other instrument trays.
- Base and insert completely nest for ease of storage.

Transpak™

NO.	DESCRIPTION	INSIDE DIMENSIONS (mm)		
		Width	Length	Depth
8600	Transpak w/Base, Lid & Insert	363	635	152
8600-00	Base & Lid			
8600-01	Base Only			
8600-02	Lid Only			
8600-03	Insert Only			

WA Distributor:

MEDICAL

■ PHONE: 1800 SONITA(766 482)

■ MOBILE: 0423 178 332

■ FAX: 08 9446 2637

SONITA MEDIC

SURGICAL

■ EMAIL: customers@sonitamedic.com.au

■ WEB: www.sonitamedic.com.au

SUPPLIES

■ OFFICE ADDRESS: Suite 76, Ground Floor 50 St. Georges Terrace, Perth WA 6000